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## ENGLISH

### 1. FUNCTION AND INDICATION

The Bleed STP+ Gel should be used in procedures where the acceleration process of the natural hemostasis of blood is required. The product is indicated for surgical procedures.

#### ⚠ WARNING

This product is of restricted use to the surgeons dentists/doctors and professionals of the rescue area.

#### ⚠ OBSERVATION

The product is absorbed by the organism in a period between 24 - 48 days.

### 2. COMPOSITION

Natural plant polysaccharide, carboxymethylcellulose and water for injectables.

### 3. ACTION MECHANISM

The formulation of the Bleed STP+ Gel product was especially designed for its use in procedures where the control of hemostasis is required. When in contact with blood, the vegetable polysaccharide absorbs the most fluid part of this (plasma) to form a gel. The extremities of this gel concentrate platelets, red blood cells and other blood constituents, accelerating the natural process of hemostasis.

### 4. PRESENTATION FORMS

Bleed STP+ Gel is sold in the following Kits:

#### • Kit with 3 g:

- 1 Applicator (syringe type) with 0.40 g of Bleed STP+ Gel (powder phase);
- 1 Applicator (syringe type) with 3.50 g of Bleed STP+ Gel (liquid phase);
- 1 Syringe adaptor;
- 1 Cannula;
- 1 User manual (available electronically).

#### • Kit with 9 g:

- 3 Applicators (syringe type) with 0.40 g of Bleed STP+ Gel (powder phase) each;
- 3 Applicators (syringe type) with 3.50 g of Bleed STP+ Gel (liquid phase) each;
- 3 Syringes adaptor;
- 3 Cannulas;
- 1 User manual (available electronically).

### 5. STORAGE AND TRANSPORT

Keep away from direct light and heat, and away from the reach of children. Store product at room temperature.

### 6. STERILE PRODUCT


The product is supplied sterile, therefore packs should be sealed before use. Bleed STP+ Gel is sterilized with gamma rays.

#### ⚠ WARNING

The material cannot be re-sterilized. In case that it has material surplus after the use, this will have to be discarded.

### 7. USE INSTRUCTIONS

1. Open outer packaging (quaternary) and check whether tertiary (blister) is intact;
2. Open tertiary packaging (blister) when preparing surgical table;
3. Open secondary packaging (surgical grade paper) only at time of product application;
4. Remove powder-phase and liquid-phase syringes covers;
5. Couple the two syringes with the aid of the syringe adaptor (as shown below);

6. Start the to-and-fro system, start by adding all the liquid to the powder. This movement should be repeated 20 times;
7. The final product, that is, the combination of powder and liquid phases should be used within a maximum of two hours after mixing, and it should not be used on another patient or surgery.
8. Apply product directly on the site where bleeding occurs, takes place and, if necessary, apply slight pressure with a gauze. If bleeding persists, remove excess of product and repeat application. After 30 seconds, in case gauze is used, wet it either with water for injectables or a physiological saline solution and removing it, observe whether bleeding has stopped.

### 8. GENERAL SAFETY INSTRUCTIONS

- Use is restricted to Dental Surgeons/ Surgeons and rescue professionals;
- Product application should be restricted to the treatment area;
- This product should be handled by an operator wearing protective gloves;
- Once hemostasis is obtained, the excess of Bleed STP+ Gel should be removed from the application site by irrigation and aspiration, especially when used in bone foramens, bone regions, spinal cord and/or in the nerve or optic chiasm. Product increases in volume when in contact with liquids. In case product is used as a hemostatic agent, the Bleed STP+ Gel remaining after the application should be removed of dry material in excess reduces the chances of product interfering with normal physiological functions and/or causing necrosis by the compression of adjacent tissues.
- The safety and efficacy of Bleed STP+ Gel were not clinically tested on children or pregnant women. It has been observed that amylase activity is much lower in newborns and babies up to the first ten months of life; therefore, Bleed STP+ Gel may take more than 40 days to be absorbed in such population;
- Bleed STP+ Gel should be used with caution in the case of infections, or if used in contaminated areas of the body. In case of any sign of abscess or infection after product application, a new drainage intervention may become necessary;
- Bleed STP+ Gel should not be used to control postpartum hemorrhage or menorragia.

### 9. SIDE EFFECTS

There are no side effects reports for this product.

### 10. CONTRAINDICATIONS

- Do not mix this substance with vasoconstrictors: epinephrine, racemic epinephrine, iron chlorides, ferrous sulfate, aluminum chloride, etc;
- Do not inject or use Bleed STP+ Gel in blood vessels as embolization or death may occur;
- Bleed STP+ Gel is not recommended as the primary treatment of coagulation disorders;
- Tests were not conducted by using Bleed STP+ Gel on bone surfaces before prosthetic materials placement with adhesives; therefore, such use is not recommended;
- Do not apply more than 100 g of Bleed STP+ Gel on diabetic patients as, according to calculations, quantities greater than this may influence glucose levels.

### 11. PRECAUTIONS

- Use Bleed STP+ Gel until the expiry date;
- Bleed STP+ Gel was designed to be used in the gel state, and its hemostatic efficacy is reduced when in contact with saline solutions or antibiotics;
- Bleed STP+ Gel should not be left in the renal pelvis or in ureters in urological procedures, in order to eliminate possible kidney stone formation focuses.

### 12. LIMITATIONS

In case this product is used in cranial operations, remove as much of the product as possible before final suture.

### 13. EXPIRY

This product is valid for 2 years from date of manufacture. After opening the surgical grade paper packaging (tertiary and secondary), it must be used immediately.

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EN

Pursuing practicality and easy access to the User Manual of DMC® products, in accordance with IN nº 4/2012 established by ANVISA, the document is available through the electronic address [www.dmcgroup.com.br](http://www.dmcgroup.com.br), according to the download instructions available in it. Important: Upon receiving the file, check that the User Manual revision matches the one indicated on the package.

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